


Experimental Design *Quality Control, Normalization*

Agnes Paquet
SincellTE, Roscoff 05/02/2019



Single RNAseq workflow: bioinformatics point of view

- What technique should we use to generate the data ?
 - Plate based / droplets
 - Full length / 3' counting with UMI
 - **UNDERSTAND THE BIAS**
- Experimental design
 - Sequencing strategy
 - UMI design
 - Spike-ins
 - Sequencing strategy?
 - Number of cells
 - Samples: Practical considerations
 - Types /number of samples
 - Cell preparation -> *confounding*
 - Budget



Technology Overview

Table 2 The advances of single-cell capture methods

Methods	Advantage	Drawback	Application
Mouth pipetting	Low cost	Time consuming	Rare sample
Laser capture microdissection	Visualization	Time consuming	Specific target
Flow cytometry	Marker selection	Require sorting	MARS-seq
Microwell platform	High throughput	mRNA capture rate	Cyto-seq
Microdroplet platform	High throughput	mRNA capture rate	Drop-seq, inDrop
Fluidigm C1 platform	Automatic library prep	High cost	qPCR, mRNA-seq
DEPArray	Visualization	High cost	Specific target

Ye et al, Journal of Hema and Onco 2017

Table 1 Brief overview of scRNA-seq approaches


Protocol example	C1 (SMARTer)	Smart-seq2	MATQ-seq	MARS-seq	CEL-seq	Drop-seq	inDrop	Chromium	SEQ-well	SPLIT-seq
Transcript data	Full length	Full length	Full length	3'-end counting	3'-end counting	3'-end counting	3'-end counting	3'-end counting	3'-end counting	3'-end counting
Platform	Microfluidics	Plate-based	Plate-based	Plate-based	Plate-based	Droplet	Droplet	Droplet	Nanowell array	Plate-based
Throughput (number of cells)	10^2-10^3	10^2-10^3	10^2-10^3	10^2-10^3	10^2-10^3	10^3-10^4	10^3-10^4	10^3-10^4	10^3-10^4	10^3-10^5
Typical read depth (per cell)	10^6	10^6	10^6	10^4-10^5	10^4-10^5	10^4-10^5	10^4-10^5	10^4-10^5	10^4-10^5	10^4
Reaction volume	Nanoliter	Microliter	Microliter	Microliter	Nanoliter	Nanoliter	Nanoliter	Nanoliter	Nanoliter	Microliter
Reference	[63]	[57]	[39]	[10]	[64]	[45]	[46]	[47]	[101]	[38]

Haque et al, Genome Medicine 2017

- Throughput: 10s vs. 1000s of cells?
- Full length protocols: required for splicing, inferring CNV, TCR/BCR profiling

- UMI design
- Use of Spike-ins
- Discuss about sequencing design
 - Number of cells
 - Sequencing depth


Unique Molecular Identifier (Islam et al., Nature Methods, 2014)



UMIs : Kivioja, T. et al. Counting absolute numbers of molecules using unique molecular identifiers. Nat Meth 9, 72-74 (2012)
UMIs for single cell transcriptome: Islam, S. et al. Quantitative single-cell RNA-seq with UMI . Nat Methods 11, (2014).

UMI design and bias correction


- UMI-based protocols allow for PCR bias correction
- Improved accuracy of gene expression measures




Chen, Genome Bio 2018

UMI Design Bias

- Design limits
 - N=4-10bp barcodes -> 4^N possible UMIs
 - N=5 -> 1024 UMIs available
 - N=10 -> 1,048,576 UMIs available
- GC content bias




Svensson V, et al 2017



Arguel MJ et al, 2017

Handling of sequencing errors






- Same UMI does not necessarily mean same molecule
 - Biases in UMI frequency and short UMIs
 - Correction for UMI saturation:
 - e.g. Grün, 2014
 - Different UMI does not necessarily mean different molecule
 - Sequencing errors
 - Different transcript does not necessarily mean different molecule
 - Mapping errors/multi-mapping
- Error correction using edit distance (ed=1 standard for 8-10bp UMI)
- Ref: UMI-tools, Smith T, Genome Res 2017

Spike-ins

- Spike-ins are molecules that are added in known concentration to the library
- Used to assess protocol accuracy and reproducibility
- ERCC
 - 92 bacterial RNA species, different lengths, GC contents
 - 22 abundance levels, 2 mixes for fold-change accuracy assessment
- SIRV
 - 69 artificial transcripts
 - Mimic human genes
 - Main difference: Used for isoforms detection


Spike-ins use in scRNA-seq

- Estimate protocol capture efficiency
 - How many of the spiked molecules did we detect?
- Comparison of protocols performance
 - Level of detection in low expressed genes
 - See Svensson V. et al, 2017
- Estimate technological noise
 - Help for detection of highly variable genes
- Issue 1: spike-ins behave differently than endogenous genes
 - May introduce more bias
- Issue 2: Spike-ins can't be used in droplet assays
 - Even incorporation in all droplets
 - Reads will be used to sequence only spike-ins



Experimental Design


- We have a question
- We have selected a protocol
- How many samples?
- How many cells?
- How many reads/cell?
- How do we combine all this to minimize batch effect?



Estimating the required number of cells / sequencing depth

- Number of cells required
 - Do we have a lot of cells to begin with?
 - Are we looking for rare cells (probability estimation)?
- WARNING: doublet rate increases with higher cell numbers in droplet assays.
- Sequencing depth
 - What are the limits of my sequencer? (Novaseq or NextSeq)
 - Minimal number of reads for droplets: 50,000 reads/cells
 - Do the cells have lots of RNA ?
 - *Think about sequencing saturation*
 - *Think about dropouts generation*

Zheng 2017



Example 1: PBMC (small cells, some don't have a lot of RNA)

-Target: 5,000 cells

- 1 sample, NextSeq High 75 (~400millions reads / run)

Estimated Number of Cells

6,388

Mean Reads per Cell

61,047


Median Genes per Cell

1,343

Sequencing


Number of Reads	389,969,360
Valid Barcodes	97.9%
Reads Mapped Confidently to Transcriptome	52.5%
Reads Mapped Confidently to Exonic Regions	54.6%
Reads Mapped Confidently to Intronic Regions	21.4%
Reads Mapped Confidently to Intergenic Regions	3.8%
Reads Mapped Antisense to Gene	3.8%
Sequencing Saturation	82.5%

Cells




Estimated Number of Cells	6,388
Fraction Reads in Cells	92.3%
Mean Reads per Cell	61,047
Median Genes per Cell	1,343
Total Genes Detected	21,143
Median UMI Counts per Cell	4,480

Sequencing Saturation






Median Genes per Cell




Example 2: Nasal epithelium brushing (cells with lots of RNA)

-Target: 5,000 cells

- 2 samples, NextSeq High 75 (~400millions reads / run)




Number of cells: example of the 1.3millions cells dataset



Bhaduri A, BiorXiv 2017


Summary of technical design

- Discuss about sequencing depth with the biologist
- If the sequencing is too shallow, the statistical analysis may not be robust
 - Worst case scenario: you can't even find the biologist favorite gene
- More cells is not always better
- Sequencing depth should be the same for all samples






Single RNAseq workflow: bioinformatics point of view

- What technique should we use to generate the data ?
 - Plate based / droplets
 - Full length / 3' counting with UMI
 - UNDERSTAND THE BIAS
- Experimental design
 - Sequencing strategy
 - UMI design
 - Spike-ins
 - How to sequence
 - **Samples: Practical considerations**
 - Types /number of samples
 - Cell preparation -> *confounding*
 - Budget



What about experimental confounding factors ?


- Most scRNA-seq are performed 1 sample at a time
 - Dissociation is difficult, sample are collected 1 by 1,...
 - Technological aspects vary too (seq depth, number of cells captured)
- Several studies report evidence for strong batch effects




Hicks , 2017

Dissociation induced genes

- Van den Brick S, Nat Method 2017




Perfect study design



- Balanced design will be hard to achieve for practical reasons

- Multiplexing :
 - Natural SNPs (demuxlet)
 - Expression of Xist/ChrY
 - Cell -hashing**



Stoeckius, 2018

ARTICLE

<https://doi.org/10.1038/s41586-018-0590-4>

Single-cell transcriptomics of 20 mouse organs creates a *Tabula Muris*


The Tabula Muris Consortium*



Cell

Mapping the Mouse Cell Atlas by Microwell-Seq

Graphical Abstract



Resource

Authors

Xiaoping Han, Renying Wang,
Yincong Zhou, ..., Guo-Cheng Yuan,
Ming Chen, Guoji Guo

Correspondence


xhan@zju.edu.cn (X.H.),
ggj@zju.edu.cn (G.G.)

In Brief

Development of Microwell-seq allows construction of a mouse cell atlas at the single-cell level with a high-throughput and low-cost platform.

Mouse Cell Atlas Summary


- > 400,000 cells
- >50 mouse tissues and cultures
- > 800 cell types identified
based on 60,000 good QC cells



ISSUE: how do you deal with >100,000 cells?


Tabula Muris Summary

- Over 100,000 cells
- 20 organs
- Double design:
 - Shallow profiling using droplets
 - FACS + full length profiling



MCA Lung data (6940 cells)


Han et Al, Cell (2018)



MCA Lung data (6940 cells)

Han et Al, Cell (2018)


A



Gene expression and cell type markers available on : <http://bis.zju.edu.cn/MCA/gallery.html?tissue=Lung>

Single-cell transcriptomics of 20 mouse organs creates a *Tabula Muris*

The Tabula Muris Consortium*



SMART-SEQ + FACS


Lung	Trachea
1620 cells	1392 cells

10X Microfluidic droplet

Lung	Trachea
5449 cells	11269 cells


TM Lung 10X data (5449 cells)

QC metrics



TM Lung 10X data (5449 cells)


T-SNE of cell types



➤ 15 cell types (8 immune)


- B cell n = 205
- ciliated cell n = 41
- Clara cell n = 5
- dendritic cell n = 225
- endothelial cell n = 425
- leukocyte n = 151
- macrophage n = 456
- mast cell n = 22
- mesothelial cell n = 24
- monocyte n = 145
- natural killer cell n = 832
- stromal cell n = 2534
- T cell n = 246
- type II pneumocyte n = 89
- unknown n = 49

TM Lung 10X data (5449 cells)




TM Lung SMART-Seq data (1620 cells)

QC metrics



TM Lung SMART-Seq data (1620 cells)

T-SNE of cell types





➤ **16 cell types (7 immune)**

- B cell **n = 55**
- ciliated cell **n = 14**
- Clara cell **n = 13**
- dendritic cell **n = 69**
- endothelial cell **n = 738**
- epithelial cell **n = 9**
- leukocyte **n = 31**
- lung neuroendocrine cell **n = 2**
- macrophage **n = 69**
- mesothelial cell **n = 2**
- monocyte **n = 65**
- natural killer cell **n = 36**
- stromal cell **n = 366**
- T cell **n = 55**
- type I pneumocyte **n = 2**
- type II pneumocyte **n = 94**

TM Lung SMART-Seq data (1620 cells)

T-SNE of batches






MCA (query)
on
TM (ref)

35%

90%

TM (query)
on
MCA (ref)



Comparing MCA and TM Lung
10X datasets

Conclusion

	MCA lung 10X	TM lung 10X	TM lung SMART-Seq
Nb of Cells	6940	5449	1620
Nb of cell types	30	15	16
Sequencing depth (mean of detected genes)	764	1200	2000
% of Dropouts	96%	93%	89%

Shared cell types

12

12

Well mapped cells

High depth on Low depth

90%

80%


Well mapped cells

Low depth on High depth

35%

45%

Comparison of the Mouse Atlases




Tabula Muris, 2018

References

- Svensson V et al, Power analysis of single-cell RNA-sequencing experiments, Nature Methods 2017
- Baran-Gale et al, Experimental design for single-cell RNA sequencing, Brief Functional Genomics 2017
- Tung PY et al, Batch effects and the effective design of single-cell gene expression studies, Science Reports 2017
- Arguel MJ et al, A cost effective 5 selective single cell transcriptome profiling approach with improved UMI design, Nuc Acid Res, 2017
- Chen et al, UMI-count modeling and differential expression analysis for single-cell RNA sequencing, Genome Biol 2018
- Grün D et al, Validation of noise models for single-cell transcriptomics, Nat Method 2014
- Ziegenhain C et al, Comparative Analysis of Single-Cell RNA Sequencing Methods, Molecular Cell 2017
- Hicks SC, Missing data and technical variability in single-cell RNA-sequencing experiments; Biostatistics 2017
- Kang HM et al, Multiplexed droplet single-cell RNA-sequencing using natural genetic variation, Nature Biotech 2017
- Stoeckius M, Cell 'hashing' with barcoded antibodies enables multiplexing and doublet detection for single cell genomics, BiorXiv 2017
- Van den Brick S, Single cell sequencing reveals dissociation-induced gene expression in tissue subpopulations, Nat Method 2017



Standard Data Analysis Pipeline



Todorov, 2018


Preprocessing Steps

- Filtering of poor quality cells
 - Number of genes/UMI detected
 - % mitochondrial genes...
- Remove doublets
 - doubletFinder, scrublet
- Check for background issues
 - SoupX




Young MD, BiorXiv 2018

BEFORE SoupX



AFTER SoupX




Warning: the software requires manual tuning.

NORMALIZATION

- Process of **identifying** and **removing** systematic variation not due to real differences between RNA treatments i.e. differential gene expression.

- Cell-specific effects

- Gene-specific effects



C

	Cell-specific effects	Gene-specific effects	Not removed by UMIs
Sequencing depth	✓		✓
Amplification	✓	✓	
Capture and RT efficiency	✓	✓	✓
Gene length		✓	
GC content	✓	✓	✓
mRNA content	✓		✓

Vallejos CA, 2017

- Gene-specific effects
 - within cell: GC content, gene length
- Cell specific effects
 - Aim: make count distributions comparable
- Sample/Technology-specific effects -> Data Integration
 - Batch effects (BAD)
 - Between samples variability (GOOD)


Bulk RNAseq normalization

- RPKM/FPKM/TPM/CPM (Reads/Fragments per kilobase of transcript per million reads of library)
 - Normalize for sequencing depth and transcript length at the same time
 - > ok if you have full length data
- Global scaling
 - Eg. Upper Quartile
 - If we have too many zeros, the SF will be off
- Size factors calculation
 - Estimation of library sampling depth
 - DESeq2, edgeR TMM
 - Suppose that **50%** of genes are not DE
 - If we have too many zeros, the SF will be off
- These methods don't work well for single-cell data
 - TPM/CPM can be bias by a small number of genes carrying most of the signal
 - Quantile based methods are limited: large number of zeros -> scale factor = 0

scRNA-seq: 3 levels of normalization

- Gene-specific effects
 - within cell: GC content, gene length
 - ***Not really accounted for in droplet assays***
- Cell specific effects
 - Aim: make count distribution comparable
 1. Global scaling
 2. scRNA-seq specific method from scater/scrn package
 3. Others
- Sample/Technology-specific effects -> Data Integration
 - Batch effects (BAD)
 - Between samples variability (GOOD)


- Hypotheses:
 - Cell populations are homogenous
 - The RNA level is similar in all cells
- Choice of the scaling factors
 - Median UMI counts
 - 10,000 default in Seurat / Cell Ranger




- In practice
 - Hypotheses are not always verified, but lots of people use this method anyway

Estimation of size factors using deconvolution

- Alternative method to compute the size factors
- Pool cells to reduce the number of zeros
- Estimate the size factors for the pool
- Repeat many time and use deconvolution to estimate each cell size factor
- Implemented in **scater**/**scran** packages



Lun, 2016



Vallejos C, 2017


Other methods are available

- Normalization included in the statistical model
 - SCDE, Monocle, MAST,...
- Normalization based on spike-ins or invariant genes
 - BASiCs
- Can we be more creative?

Common Approach:

Normalizing independent of cell types

Observed Count Matrix
Cells



Normalization

To mean/median library size
Downsampling
BASiCs with spike-ins/ERCCs



Clustering Cells



Downstream
Analysis

Problems:


- Dropouts not resolved
Zeros remain zero!
- Removes biological stochasticity specific to cell type
- Leads to improper clustering; Biased downstream analysis

Azizi, 2017

Integration of imputation, pre-processing and clustering

- SCONE (R package)

Cole M, Risso D (2018). scone: Single Cell Overview of Normalized Expression data. R package version 1.4.0.



Nir Yosef

- BISCUIT (R package)

Bayesian Inference for Single-cell Clustering and Imputing
Elham Azizi, 2017

Normalization for other biological factors


- Known or unknown variation
 - Cell cycle, number of genes detected, % mitochondrial genes...
- Regression methods provided to account for know factors
 - Seurat
- Latent variable models to estimate and remove unknown bias
 - scLVM

scRNA-seq: 3 levels of normalization

- Gene-specific effects
 - within cell: GC content, gene length
- Cell specific effects
 - Aim: make count distribution comparable
 1. Global scaling
 2. scRNA-seq specific method from scater/scrn package
 3. Others
- Sample/Technology-specific effects -> Data Integration
 - Batch effects (BAD)
 - Between samples variability (GOOD)

Why do we need data integration methods?

- In practice: single cell techniques are biased
 - Variations between samples can be huge
 - donor effect +/- sampling effect
 - Samples may be processed using different technologies
- Combining datasets and applying cell-level normalization might not be enough to remove this bias





Misharin, BiorXiv 2018

For differential analysis:

- > Choose a framework where you can add a batch term in your statistical model (e.g.: MAST, DESeq2, limma,...)


For other analyses:

- We need a method that will “merge” our datasets and remove the unwanted variation
- Non-linear transformation of cells in different proportions
- Aligns datasets from different technologies and species



Current Data Integration Method

- MNN: Haghverdi, 2018
- Harmony, Korsunsky BiorXiv 2018
- Seurat V3, Stuart BiorXiv 2018
- ComBat (sva)
- No gold-standard yet
- Performance assessment?
 - Visual inspection
 - kBET (Büttner, 2019)
 - Other metrics?



Haghverdi, 2018

References

- Vallejos CA, Normalizing single-cell RNAsequencing data: challenges and opportunities, Nat Method 2017
- **Scater**: Lun A, Pooling across cells to normalize single-cell RNA sequencing data with many zero counts, Genome Biology 2016
- **Seurat**: Butler et al. Integrating single-cell transcriptomic data across different conditions, technologies, and species. Nature Biotechnology (2018).
<https://satijalab.org/>
- https://bioconductor.org/help/course-materials/2017/BioC2017/Day2/InvitedSpeakers/Biscuit_Azizi.pdf
- Haghverdi, L., Batch effects in single-cell RNA-sequencing data are corrected by matching mutual nearest neighbors. Nat. Biotechnol. 2018.
- kBET: Maren Büttner et al, A test metric for assessing single-cell RNA-seq batch correction, Nat Methods 2019

Comparing datasets with SC-Map

Kiselev et al., Nature Methods (2018)

SC-Map = R package with a label-centric approach, focused on trying to identify equivalent cell-types across datasets by comparing individual cells or groups of cells.

Method used = sc-map cluster

- 3 steps :**
- 1) Selection of most informative genes (cell types markers)
 - 2) Compute expression median for selected genes in all cells of each cluster
 - 3) Correlation tests between each query cells and reference expression profiles

