



Single-cell RNA-seq analysis pipelines using the Eoulsan Workflow engine


**Laurent Jourden, Nathalie Lehmann,
Morgane Thomas-Chollier**

*Computational systems biology - IBENS
& Genomic Platform Paris Centre
(Stéphane Le Crom)*

Question : how to process a scRNA-seq dataset ?



Question : how to process a scRNA-seq dataset ?



Eoulsan : a workflow engine designed for bulk RNA-seq

- Build on existing software




Eoulsan

Jourdren L, Bernard M, Dillies MA, Le Crom S.
Bioinformatics. 2012

- Modularity and reproducibility ( containers)




Users (bioinformatician)




- Input/output complexity is hidden
- User just need to focus on the workflow design

Users (bioinformatician)


Developer



Design file



Configuration file



Workflow file
(XML)



FASTQ



Program 1 (parameters)

Program 2 (parameters)

Program 3 (parameters)




Expression Matrix

- Input/output compatibility
- Ensures each program is encapsulated

- Input/output complexity is hidden
- User just need to focus on the workflow design



Result : adaptation of Eoulсан for scRNA-seq

- Identification and evaluation of state-of-the-art tools to integrate as modules to Eoulсан



Result: original development of a filter on cell quality


- Problem : detect and remove cells of low quality (=noise in data, artefactual cell variability)
 - Evaluate proportion of reads mapping to mitochondrial genes (<10%)
 - Threshold on mapped reads (>500,000/cell) and detected genes (>1500)
 - Filter on sequencing **saturation plots**



10x Genomics uses Unique Molecular Identifiers for 3' end transcriptomes

- 1 barcode per cell (16bp)
- 1 UMI per transcript (10bp random barcodes added to transcripts)


Final Library Construct



10x Genomics

- In practice : don't get the full transcript (3' end only)

smart-seq2



Read quality filtering	FASTQC
Alignment (mapping)	STAR
Alignment filtering	Picard, Samtools
Expression estimation	HT-Seq count
Expression merger	Original program
Cell quality filtering	Original program
Normalization across cells	Scran

Expression level

gene 1
gene 2
gene 3
.
.

TSV / SingleCellExperiment (R object) formats

What should be changed ?


Barcode + UMI



Result : adapt Eoulsan scRNA-seq to support



Nathalie Lehmann
Engineer CNRS



Read quality filtering

FASTQC

Cell identification/filtering

UMI-tools

UMI processing

UMI-tools

Alignment (mapping)

STAR

Alignment filtering

Samtools

Expression estimation

HT-Seq count, featureCounts, UMI-tools count

Expression merger

Original program


Cleaning matrix

Original program

Normalization across cells


Scran

Expression level




Application: Eoulsan scRNA-seq on 10x GENOMICS® dataset

- Dataset from Vassili Soumelis (Institut Curie, Paris)
- Human immune cell
(plasmacytoid dendritic cell)
- 7 samples
- ~1.000 cells/sample



Eoulsan scRNA-seq is faster than Cell ranger



- Lower memory usage with Eoulsan

Euclsan scRNA-seq also supports some downstream analyses

Reads



Expression
matrix

- Biological-driven analyses
- Still very active methodological developments
- Alternative options

Downstream analyses



Differential
Expression

SCDE

Clustering

Seurat

Lineage
analysis

Monocle2

- **Eoulsan**

<https://github.com/GenomicParisCentre/eoulsan>



- Install
- New user's tutorial (bulk, smart-seq, 10x)
- Complete documentation