

# Singel Cell RNA-seq: Technologies and Experimental Approaches

Kévin Lebrigand  
UCAGenomix, Nice-Sophia-Antipolis

✉ [lebrigand@ipmc.cnrs.fr](mailto:lebrigand@ipmc.cnrs.fr)



🐦 @kevinlebrigand

Roscoff, 18 juin 2018





# Methods of the year 2013 (Nature Methods)

- 2007: Next Generation Sequencing
- 2008: Super-resolution microscopy
- 2009: Induced pluripotency
- 2010: Optogenetics
- 2011: Genome editing with engineered nucleases (TALEN)
- 2012: Targeted proteomics
- 2013: Single cell sequencing
- 2014: Light-sheet fluorescence microscopy
- 2015: Cryo-EM or electron cryomicroscopy
- 2016: Epitranscriptome analysis




# Why single cell profiling?

Stop measuring gene average as in bulk,  
Population sequencing yields average values,  
Changes in subpopulation might remain undetected,




# Why single cell profiling?


## 1. Assess cell-to-cell heterogeneity




## 2. Map cell trajectories




## 3. Dissect transcriptional mechanics



## 4. Infer gene regulatory networks




Network inference



Liu et al., F1000, 2016




# Why single cell profiling?



- Understanding heterogeneous tissues,
- Identification of **cell types** and **cell states**
- Identification and analysis of **rare cell types**
- **Changes in cellular composition**
- Dissection of temporal changes

Sandberg et al.,  
Nature Methods, 2014

# Single cell analysis: the context




- Highly efficient library preparation techniques
- Elimination of PCR amplification bias
- Spike-in ERCC molecules to evaluate yield and capture efficiency
- Use of Unique Molecular Identifiers to monitor the number of molecules:
  - *Random multimers to reduce cloning biases during small RNAs profiling (Jayaprakash, NAR, 2011)*
  - *Detection of rare DNA mutations (Ståhlberg et al, NAR, 2016)*
  - *Improved accuracy of molecule counting*


# ERCC spike-ins

ERCC (Externals RNA Controls Consortium) spikes is 92 polyadenylated RNA molecules

- different sequences and lengths,
- no homology with known genome sequences,
- 2 mixes (mix1 and mix2) with the 92 sequences but in different amounts
- relative amounts of spikes in mix1 and mix2 from 1 to 10e6




Scatter plot of read counts for each ERCC transcript in 2 RNA-seq libraries  
Jiang et al., Gen.Res. (2015)



Single Cell capture efficiency  
Arguel et al., NAR (2016)


# Unique Molecular Identifier (Islam et al., Nature Methods, 2014)



UMIs : Kivioja, T. et al. Counting absolute numbers of molecules using unique molecular identifiers. Nat Meth 9, 72-74 (2012)  
UMIs for single cell transcriptome: Islam, S. et al. Quantitative single-cell RNA-seq with UMI . Nat Methods 11, (2014).




# Unique Molecular Identifier (Islam et al., Nature Methods, 2014)



UMIs : Kivioja, T. et al. Counting absolute numbers of molecules using unique molecular identifiers. Nat Meth 9, 72-74 (2012)  
UMIs for single cell transcriptome: Islam, S. et al. Quantitative single-cell RNA-seq with UMI . Nat Methods 11, (2014).


# Unique Molecular Identifier (Islam et al., Nature Methods, 2014)




UMIs : Kivioja, T. et al. Counting absolute numbers of molecules using unique molecular identifiers. Nat Meth 9, 72-74 (2012)  
UMIs for single cell transcriptome: Islam, S. et al. Quantitative single-cell RNA-seq with UMI . Nat Methods 11, (2014).

# To measure sequences in individual cells, need methods that capture one cell at a time


**MICROPIPETTING  
MICROMANIPULATION**




**LASER CAPTURE  
MICRODISSECTION**




**FACS**



**MICRODROPLETS**




**MICROFLUIDICS  
e.g. FLUIDIGM C1**




# Manual Cell sorting





## Manual dissociated cells



fluorescence  
microscope




~100 cells




[http://www.nature.com/nprot/journal/v6/n5/images\\_article/nprot.2011.322-F2.jpg](http://www.nature.com/nprot/journal/v6/n5/images_article/nprot.2011.322-F2.jpg)







## LCM: laser capture microdissection



<http://www.genomemedicine.com/content/figures/gm247-2-l.jpg>




## FACS: fluorescence activated cell sorting




[http://www.flowlab-childrens-harvard.com/yahoo\\_site\\_admin/assets/images/principle123.285181420\\_std.gif](http://www.flowlab-childrens-harvard.com/yahoo_site_admin/assets/images/principle123.285181420_std.gif)

# Microfluidics: Fluidigm C1




- Limiting factor is size of capture chambers (96 chips: 5-10 $\mu$ m, 10-17 $\mu$ m, >17 $\mu$ m)
- 800 cells chip (10-17  $\mu$ m diameter cells)
- ScriptHub: protocols for running SMARTer, SmartSeq2, CEL-seq, STRT, ...
- openApp Chip for custom protocols development






- water-in-oil droplets to capture cells with pre-barcoded (BC+UMIs) beads,
- really really high throughput (40,000 cells),
- No cell size issue





# Scaling Single Cell Transcriptomics (Svensson et al., Nature Methods, 2017)




## Protocol

- Total RNA is isolated and fragmented,
- Converted to cDNA by using an **oligodT** primer with a specific anchor sequence,
- **Second strand synthesis using a polyT primer** with another anchor sequence,
- **PCR amplified** from primers against the two anchor sequences.

## Drawback


- Premature termination of RT reduces transcript coverage at the 5' end
- Introduction of a polyA tail in addition to its own polyA sequence at the **3' end** of the input RNA causes a **loss of strand information** in the resulting double-stranded cDNA



# STRT: Single cell tagged reverse transcription, Islam, Nat.Prot. (2012) @ Linnarsson's lab

## Protocol

- based on **template switching**,
- **5' end cDNA tagged N5 UMI**,
- biotin is introduced at both the 3' and 5' ends via the use of biotinylated primers.
- enzymatic cleavage leads to the selection of only the 5' fragments for library construction.
- sequencing and analysis shows **5' read bias**



# SMART-Seq, Ramskold (2012), SMART-seq2, Picelli (2014) @ Sandberg's lab


SMART= Switching Mechanism at the end of the 5'-end of the RNA Transcript

## Protocol

- Based on template switching mechanism,
- Anchor a 5' universal seq. along with Locked nucleic acid by reverse transcription,
- cDNA is then **PCR amplified**,
- Tagmentation is used to construct libraries,
- Generate **full transcript coverage**

## Drawback

- No UMIs






# Cell expression by linear amplification and sequencing (CEL-seq, CEL-seq2)


## Protocol

- OligodT primer containing the 5' Illumina adaptor, a cell barcode, and a T7promoter (**CEL-seq2 add UMI**),
- RT and second-strand synthesis,
- cDNA from all the cells is pooled and amplified by **in vitro transcription** from the T7 promoter,
- RNA fragmentation, Illumina adaptor ligated at 3' end
- RNA is reverse transcribed, library is prepared then sequencing,
- Sequencing of the **3'terminal fragments**



# A cost effective 5' selective single cell transcriptome profiling approach

- UMIs inside Fluidigm C1, PCR bias removal (molecules counting)
- 60k transcripts in HEK (homogeneity of chemistry)
- Illumina and Ion Torrent sequencing (no paired-end required)
- Capture efficiency **26%** (ERCC spike-ins)
- Sequences 5' end of transcripts → TSS identification



Now available on ScriptHub (Fluidigm)



Single-cell mRNA Seq with Integrated Barcoding  
mRNA Sequencing

Arguel et al.

Pascal Barbry - UCA Genomix - IPMC, CNRS - University of Cote d'Azur

*Arguel et al., Nucleic Acid Research, 2016*


# Drop-seq / inDROP: droplets-based scRNA-seq (*Mascosko et al., Klein et al., CELL, may 2015*)

## Highly Parallel Genome-wide Expression Profiling of Individual Cells Using Nanoliter Droplets



Evan Z. Macosko<sup>1, 2, 3</sup>, Anindita Basu<sup>4, 5</sup>, Rahul Satija<sup>4, 6, 7</sup>, James Nemesh<sup>1, 2, 3</sup>, Karthik Shekhar<sup>4</sup>, Melissa Goldman<sup>1, 2</sup>, Itay Tirosh<sup>4</sup>, Allison R. Bialas<sup>8</sup>, Nolan Kamitaki<sup>1, 2, 3</sup>, Emily M. Mardersteck<sup>9</sup>, John J. Trombetta<sup>4</sup>, David A. Weitz<sup>5, 10</sup>, Joshua R. Sanes<sup>9</sup>, Alex K. Shalek<sup>4, 11, 12</sup>, Aviv Regev<sup>4, 13, 14</sup>, Steven A. McCarroll<sup>1, 2, 3</sup>

- droplets encapsulation of cells and barcoded beads
- 3' selective single cell RNA-seq
- 12bp cell barcode and 8bp UMI
- capture efficiency **12.5%**
- 44,808 mouse retinal cells
- identification of 39 different cell types

### Barcoded primer bead




### Drop-seq single cell analysis



1000s of DNA-barcoded single-cell transcriptomes

# 10X Genomics Chromium





# Experimental approaches comparison

Lab (Ref) Name of approach	Region sequenced	UMI		Capture efficiency (ERCC)	Amplification	Barcoding	Device	Sequencer
		Design	Counting / error correction					
This approach	5'	N4H4	Unique (a) / edit distance (b)	26%	PCR	PCR pre-fragmentation	C1	Illumina
Linnarsson (10) STRT-seq	5'	N5	Unique start (c) / percentile (d)	48% (e)	PCR	Tagmentation <i>custom barcoded transposons</i>	C1	Illumina
Linnarsson (7) STRT-seq			Unique start (c) / percentile + (df)	22% (b)				
Amit (13) Mars-seq	Close to 3'	N4	edit distance (b)	2 – 3%	Isothermal	RT	Plates	Illumina
McCarol (16) Drop-seq	Close to 3'	N8	edit distance (b)	12.8%	PCR	RT	Microdroplets	Illumina
Kirschner (17) InDrop	Close to 3'	N6	(g)	7.1%	Isothermal	RT	Microdroplets	Illumina
Yanai (15) Cel-seq	Close to 3'	No	none	6%	Isothermal	RT	Plates	Illumina
Yanai (15) Cel-seq2	Close to 3'	N6	none	22% in C1	Isothermal	RT	Plates, C1	Illumina
Smartseq2	Internal (a)	No	–	n.d.	PCR	–	Tubes	Illumina

Arguel et al., NAR, 2016